

CATALOG & PRICING

Cat. #	Description	Size
707-50	Use with HRP in chemiluminescent ELISA	50 mL
707-100		100 mL
707-200		200 mL

BACKGROUND

Enhanced chemiluminescence (ECL) substrates are used with horseradish peroxidase (HRP)-conjugated secondary antibodies to detect protein abundance in immunoassays. HRP catalyzes luminol oxidation in the presence of hydrogen peroxide, producing light that can be detected by an X-ray film, digital imager or chemiluminescence plate reader. ShQ™ ECL Substrate Kit is proprietary formulated for chemiluminescent ELISA test. Compared to ultrasensitive PiQ™ ECL, ShQ™ ECL Substrate Kit delivers superior performance by generating an extremely low background.

PRODUCT FEATURES

- Superior performance with an ultra-low background
- Suitable for detecting highly abundance analytes
- Ready-to-use kit with easy-to-follow instructions

STORAGE

Stable at 4°C for at least 1 year.

KIT COMPONENT

Equal volume of Reagents A and B

EXPERIMENTAL PROTOCOL

1. Let the ECL reagents warm up at room temperature for approximately 30 min prior to use.
2. Combine Reagents A and B at a ratio of 1:1. Vortex to mix.
3. Apply the mixture to the microplate wells.
4. Immediately read the relative light unit (RLU) of the chemiluminescence signals using a chemiluminescence plate reader.

APPLICATION TESTED

Chemiluminescent ELISA

RECOMMENDED ANTIBODY DILUTIONS

Primary antibody: 1: 5,000 - 1:10,000

Secondary antibody: 1: 10,000 - 1:20,000

VALIDATION DATA




Figure 1. Comparison of ShQ™ ECL and PiQ™ ECL sensitivity in chemiluminescent ELISA. Different concentrations of gamma-irradiated SARS-CoV-2 viruses were lysed in the Viral Lysis Buffer (Cat# VL101) and coated (150 μ L) on microplate wells. A rabbit monoclonal anti-N primary (Cat# VYN7, 1:10,000) and an HRP-conjugated highly cross-adsorbed goat anti rabbit IgG secondary antibody (Cat# 202, 1:20,000) were used to detect SARS-CoV-2 N protein. ShQ™ ECL (Cat# 707) and PiQ™ ECL (Cat# 636) Substrate Kits were used to generate chemiluminescent signals. RLU, relative light unit. Δ RLU = RLU of the well with the primary antibody - RLU of the corresponding well without the primary antibody, everything else being equal. **Note:** When the antigen is abundant, ShQ™ ECL is a better choice based on linear regression data.

Table 1. Comparison of ShQ™ ECL and PiQ™ ECL sensitivity and background noise.

	ShQ™ ECL			PiQ™ ECL		
	w/o 1st Ab	w/ 1st Ab	Δ RLU	w/o 1st Ab	w/ 1st Ab	Δ RLU
Lysis buffer	3056	2927	-129	130727	137531	6804
Lysis buffer w/o 2nd Ab	2860	2607	-253	122382	127826	5444
105 TCID50 SARS-CoV-2	2781	73695	70914	128401	8817960	8689559
52.5 TCID50 SARS-CoV-2	2987	30316	27329	134264	5748869	5614605
42 TCID50 SARS-CoV-2	2726	23527	20801	128083	4198178	4070095
21 TCID50 SARS-CoV-2	2767	10342	7575	131388	1797070	1665682
8.4 TCID50 SARS-CoV-2	2901	4143	1242	129277	635622	506345

Notes:

1. Different concentrations of gamma-irradiated SARS-CoV-2 viruses were lysed in the Viral Lysis Buffer (Cat# VL101) and coated (150 μ L) on microplate wells.
2. A rabbit monoclonal anti-N antibody (Cat# VYN7, 1:10,000) and an HRP-conjugated highly cross-adsorbed goat anti rabbit IgG secondary antibody (Cat# 202, 1:20,000) were used to detect SARS-CoV-2 N protein.
3. ShQ™ ECL (Cat# 707) and PiQ™ ECL (Cat# 636) Substrate Kits were used to generate chemiluminescent signals.
4. RLU, relative light unit. Δ RLU = RLU of the well with the primary antibody - RLU of the corresponding well without the primary antibody, everything else being equal.