

Description


Product Name	JNK antibody
Host Species	Rabbit
Clonality	Polyclonal antibody
Purification	Polyclonal antibody is produced by immunizing animals with a recombinant human JNK fusion protein. Affinity purification using immunogen.
Applications	Western blotting
Species Reactivity	Human; Mouse; Rat; Monkey; Bovine; Zebrafish
Specificity	JNK rabbit polyclonal Ab detects endogenous levels of total JNK1, JNK2 or JNK3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	JNK
Other Names	JNK2; SAPK/JNK
Concentration	1.0 mg/ml
Storage	Store at -20°C for long term preservation(recommended).Store at 4°C for short term use.

Application Details

SDS-PAGE molecular Wt: 54 kDa and 46 kDa, respectively

Western blotting: 1:500~1:2000 (recommended 1:1000)

Images



Western blot analysis of of extracts from PC12 cells, using JNK rabbit polyclonal antibody #AT1028.

