

AR antibody



Catalog No: #AT10855

Host Species: Rabbit

Clonality: Polyclonal antibody

Species Reactivity: Human; Mouse

SDS-PAGE MW: 110 kDa

Applications: Western blotting

Package Size:
#AT10855S 20 μ l
#AT10855M 50 μ l

Description

Product Name AR antibody

Purification The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity AR polyclonal antibody detects endogenous levels of total AR.

Immunogen Type Peptide

Immunogen Description Synthesized peptide derived from human AR.

Target Name AR

Other Names Androgen receptor; ANDR; DHTR

Concentration 1.0 mg/ml


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note In order to get the best results using various techniques and preparations, we recommend detecting optimal working dilutions by titration.

Application Details

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from mouse kidney tissue, using AR rabbit polyclonal antibody #AT10855.



Tel: 025-85390309; 85300038

Orders: sales@sciben.com

Support: support@sciben.com

Website: www.sciben.com