

CHOP antibody



Catalog No: #AT12987

Host Species: Mouse

Clonality: Monoclonal antibody (mAb)

Species Reactivity: Human; Mouse; Rat

SDS-PAGE MW: 27 kDa

Applications: Western blotting
Immunohistochemistry

Package Size:
#AT12987S 20 μ l
#AT12987M 50 μ l

Description

Product Name CHOP antibody

Purification Monoclonal antibody is purified by affinity isolation.

Specificity CHOP monoclonal antibody detects endogenous levels of total CHOP.

Immunogen Type
Peptide

Immunogen Description
Synthetic Peptide

Target Name CHOP

Other Names DDIT3

Concentration 1.0 mg/ml

Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Note In order to get the best results using various techniques and preparations, we recommend detecting optimal working dilutions by titration.

Application Details

Western blotting: 1:500~1:2000 (recommended 1:1000)

Immunohistochemistry (IHC): 1:50~1:200

Images



Western blot analysis of lysates from mouse liver tissues, using CHOP mouse monoclonal antibody #AT12987.



Tel: 025-85390309; 85300038

Orders: sales@sciben.com

Support: support@sciben.com

Website: www.sciben.com