

c-Myc antibody



Catalog No: #AT13292

Host Species: Mouse

Clonality: Monoclonal antibody (mAb)

SDS-PAGE MW:
57-70 kDa

Applications:
Western blotting

Package Size:
#AT13292S 20 μ l
#AT13292M 50 μ l

Description

Product Name c-Myc antibody

Purification Monoclonal antibody is purified by affinity isolation.

Specificity The c-Myc tag antibody can detect c-Myc fusion proteins.

Immunogen Type
Recombinant Protein

Immunogen Description
Recombinant Protein

Target Name c-Myc

Concentration 1.0 mg/ml


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note In order to get the best results using various techniques and preparations, we recommend detecting optimal working dilutions by titration.

Application Details

Western blotting: 1:1000~1:2000

Images



Western blot analysis of c-Myc recombinant protein, using c-Myc mouse monoclonal antibody #AT13292.



Tel: 025-85390309; 85300038

Orders: sales@sciben.com

Support: support@sciben.com

Website: www.sciben.com