

c-Myc Tag antibody



Catalog No: #AT13296

Host Species: Rabbit

Clonality: Polyclonal antibody

SDS-PAGE MW:
55 kDa

Applications:
Western blotting

Package Size:
#AT13296S 20 μ l
#AT13296M 50 μ l

Description

Product Name c-Myc Tag antibody

Purification Antibodies were purified by affinity-chromatography using epitope-specific peptide.

Specificity The antibody detects transfected proteins containing c-Myc Tag.

Immunogen Type
Peptide

Immunogen Description
Synthetic peptide

Target Name c-Myc Tag

Other Names c- Myc epitope Tag

Concentration 1.0 mg/ml


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note In order to get the best results using various techniques and preparations, we recommend detecting optimal working dilutions by titration.

Application Details

Western blotting: 1:1000~1:3000

Images



Western blot analysis of recombinant c-Myc tagged protein, using c-Myc Tag rabbit polyclonal antibody #AT13296.



Tel: 025-85390309; 85300038

Orders: sales@sciben.com

Support: support@sciben.com

Website: www.sciben.com