

Flag-Tag antibody



Catalog No: #AT15427

Host Species: Mouse

Clonality: Monoclonal antibody (mAb)

Species Reactivity: Human; Mouse; Rat

Applications:
Western blotting
Immunofluorescence
Immunoprecipitation

Package Size:
#AT15427S 20 μ l
#AT15427M 50 μ l

Description

Product Name Flag-Tag antibody

Purification Monoclonal antibody is purified by affinity isolation.

Specificity The Flag tag antibody can recognize C-terminal, internal, and N-terminal Flag-tag fusion proteins.

Immunogen Type
Recombinant Protein

Immunogen Description
Recombinant Protein

Target Name Flag-Tag

Concentration 0.5 mg/ml

Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note In order to get the best results using various techniques and preparations, we recommend detecting optimal working dilutions by titration.


Application Details

Western blotting: 1:1000~1:3000

Immunofluorescence: 1:200~1:500

Immunoprecipitation: 1:500~1:1000

Images



Western blot analysis for flag fusion protein (upper panel), and Western blot analysis of lysates from PC12 cells infected with Ad-GFP or Ad-FLAG-dn-c-Jun (bottom panel), respectively, using Flag-Tag mouse monoclonal antibody #AT15427.



Tel: 025-85390309; 85300038

Orders: sales@sciben.com

Support: support@sciben.com

Website: www.sciben.com