

GST-Tag antibody



Catalog No: #AT16241

Host Species: Mouse

Clonality: Monoclonal antibody (mAb)

Species Reactivity:
N/A

SDS-PAGE MW:
N/A kDa

Applications:
Western blotting

Package Size:
#AT16241S 20 μ l
#AT16241M 50 μ l

Description

Product Name GST-Tag antibody

Purification Monoclonal antibody is purified by affinity isolation.

Specificity GST-Tag mouse mab detects endogenous levels of C-terminal, internal, and N-terminal GST-tagged proteins.

Immunogen Type
Recombinant protein

Immunogen Description
Recombinant protein

Target Name GST-Tag

Concentration 1.0 mg/ml


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note In order to get the best results using various techniques and preparations, we recommend detecting optimal working dilutions by titration.

Application Details

Western blotting: 1:1000~1:3000 (recommended 1:2000)

Images



Western blot analysis of GST-Tag fusion protein, using GST-Tag mouse monoclonal antibody (1:2000) #AT16241.



Tel: 025-85390309; 85300038

Orders: sales@sciben.com

Support: support@sciben.com

Website: www.sciben.com