




## Antibodies for Animal Sciences

### Validated Antibodies to Assist Your Research


Researchers in agricultural and veterinary clinics have been struggling to find the appropriate antibodies for their species of interest, and we provide a range of antibodies for farm and companion animals. To guarantee the quality of our products, we test our antibodies with cell lysates of various species. Below is an example of how we validate complement c1q binding protein(C1QBP) rabbit monoclonal antibody (Cat. #1166). Please visit [our website](#) to check if we have your antibodies of interest.

#### Validation with shRNA KD



**Figure 1.** Western blotting analysis of C1QBP expression. **(A)** C1QBP expression in various cell lines. **(B)** C1QBP expression in wild-type & shRNA-knockdown HeLa cells. GAPDH served as a loading control.

#### Cross-Species Validation



**Figure 2.** Western blotting analysis of C1QBP expression in various species. Species cross-reactivity of C1QBP antibody (Cat. #1166) was tested with lysates (40  $\mu$ g) of various cell lines.

