

Description

Product Name	Phospho-p70 S6 Kinase (Thr389) antibody
Host Species	Rabbit
Clonality	Monoclonal antibody (mAb)
Purification	Monoclonal antibody is produced by immunizing animals with a synthetic phosphopeptide corresponding to residues surrounding Thr389 of human p70 S6 kinase. Affinity purification using immunogen.
Applications	Western blotting
Species Reactivity	Human; Mouse; Rat
Specificity	Phospho-p70 S6 kinase (Thr389) rabbit mAb detects endogenous levels of p70 S6 kinase only when phosphorylated at Thr389.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	p70 S6 Kinase
Other Names	KS6B1; p70-S6K; RPS6KB1; Ribosomal protein S6 kinase; polypeptide 1
Concentration	1.0 mg/ml
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

SDS-PAGE molecular Wt: 70 kDa

Recommended dilutions for Western blotting: 1:1000

Images

